

Ti-Core® Auto E

Ideal for aesthetic restorations!



Patented,* lanthanide reinforced, dual cured, resin composite core build-up material.

FEATURES

Lanthanide reinforced

Ideal composite viscosity for core build-up

Ease of preparation

Radiopaque

Dual cure formulation

Reduces chair time

Ti-Core Auto E - Vita® A2

BENEFITS

Lanthanide increases strength without the potential corrosive by-products that result from silver oxidation of silver filled ionomers. High strength for long term durability.

Enables use of old crown as a core form with minimal composite remaining on the margin. Optimizes adaptation of core material to the post.

Easy to place. Due to its similar strength to dentin, ditching and grooving of the core material is avoided.

Easily identified on x-ray.

Speeds procedure while insuring depth-of-cure. Adaptable to any technique.

Auto-mix mechanism produces consistent fusion of base and catalyst. Eliminates pad mixing and transfer. Intra oral tip allows for precise application.

For aesthetic restorations.



Cuts Like Dentin, No Ditching or Grooving

During preparation of the core, it is easy to ditch into conventional composite and glass ionomer materials due to their weaker structure. Because Ti-Core is as strong as dentin it offers resistance similar to dentin.

From the brand you trust - Ti-Core Auto E provides the strength of Ti-Core with the ease of preparation!



**EDS
Clear
Core
Forms**

* Data on file.
* U.S. Patent Nos. 5,034,433 and 5,204,398. European Patent No. EP 0 631 598 B1. Other patents pending.

Ti-Core Auto

THE ESSENTIAL CORE MATERIAL.

Ti-Core Auto E is the only lanthanide incorporated, auto mixing, composite core materials. Lanthanide reinforcing agents are bio-compatible fillers that add strength, solidity and durability. They provide a totally unique core material with properties no other composite or glass ionomer can match.

Ti-Core Auto E was developed by placing lanthanide within a Bis-GMA matrix. The results are a core material that provides a long lasting final core restoration.

TI-CORE AUTO MIX MATERIALS ELIMINATE DITCHING AND GROOVING.


During preparation of the core, it is easy to ditch into conventional composite and glass ionomer materials due to their weaker structure. Ti-Core Auto offers resistance similar to dentin. This prevents the typical ditching and grooving associated with conventional non-reinforced core materials. As a result, Ti-Core build-up material improve and simplify the preparation of the core.



Preps like dentin, no ditching or grooving



ESSENTIAL DENTAL SYSTEMS CORE FORMS



- Tab top for easy pick-up and placement
- Reference line for accurate contouring
- Tapered slides accurate sizing and easy removal
- Clear for light curing at all angles
- Polyethylene: Will not stick to core material, easy to trim.
- 4 sizes for all clinical situations

TECHNICAL FACTS:

Post Head Retention (w/ #2 Flexi-Post)
 Diametral Tensile Strength
 Fractural Strength (w/ #1 Flexi-Post)
 Radiopacity

Ti-Core Auto

880 N (±190 N) (198 lbs (±35 lbs.))
 35 MPa (±4 MPa)
 1220 N (±170 N) (274 lbs (±39 lbs.))
 Positive

ORDERING INFORMATION:

Ti-Core Auto Syringes - Dual Cured Reinforced Core Materials

Ti-Core Auto E Catalog Number 830-00
 2 auto mix syringes (each containing 4.5 gm base, 4.5 gm catalyst), 20 mixing tips, 20 intra oral tips.

Ti-Core Flow +

Ti-Core Flow + Catalog Number 832-00
 2 auto mix syringes (each containing 4.5 gm base, 4.5 gm catalyst), 20 mixing tips, 20 intra oral tips.

Auto Mix Intra-Oral Tips (ø 1.80mm OD)

20 Tips with Mixers (Fits all standard auto-mix double barrel syringes) Catalog Number 830-20

Auto Mix Intra-Oral Tips (ø 1.00mm OD)

20 Tips with Mixers (Fits all standard auto-mix double barrel syringes) Catalog Number 870-20

Essential Dental Systems Core Forms

Clear Core Forms (Box of 48, 12 ea. of 4 sizes) Catalog Number 890-00

For our complete line of patented dental products please visit our website at: www.edsdental.com