






Monarch Surface Disinfectant Wipes and Spray provide you with a fast, effective and easy way to clean and disinfect your hard surfaces throughout the office.


## Monarch Surface Disinfectant

- Fast, broad-spectrum efficacy: ethanol based, bleach free formula kills TB, viruses\*, bacteria and fungi in just 1 minute.
- Cleaning and disinfecting in one step saves time and money.
- Effective against odor-causing bacteria.
- Compatible with stainless steel, glass, porcelain, aluminum, titanium, vinyl, nylon, polypropylene, PVC, and Teflon®.
- Fragrance-free.
- Can be used on Air Techniques CamX® cameras.

A study performed by the **Dental Advisor** examined the effects caused by an accelerated 5-year interaction between *Ultraleather* and *Dauphine* dental chair upholstery and Monarch Surface Disinfectant Wipes. No visible changes in appearance were observed for either fabric treated with Monarch Surface Disinfectant Wipes.



**Figure 1a and b:** Microscopic view of a) treated and b) untreated *Ultraleather* upholstery



**Figure 2a and b:** Microscopic view of a) treated and b) untreated *Dauphine* upholstery



The ethanol-based, bleach-free formula disinfects hard, non-porous surfaces killing a broad spectrum of pathogens.

## Effective Against a Broad Spectrum of Pathogens in One Minute

Bactericidal	Virucidal*	Tuberculocidal	Fungicidal
<ul style="list-style-type: none"> <li>Acinetobacter baumannii</li> <li>Burkholderia cepacia</li> <li>Campylobacter jejuni</li> <li>Vancomycin Resistant Enterococcus faecium [VRE]</li> <li>Escherichia coli O157H7 [E coli]</li> <li>Klebsiella pneumonia</li> <li>Listeria monocytogenes</li> <li>Methicillin resistant Staphylococcus aureus [MRSA]</li> <li>Pseudomonas aeruginosa</li> <li>Salmonella enterica</li> <li>Staphylococcus aureus [Staph]</li> <li>Streptococcus pyogenes [Strep]</li> </ul>	<ul style="list-style-type: none"> <li>Hepatitis B</li> <li>Hepatitis C</li> <li>Herpes Simplex virus type 1</li> <li>Coronavirus</li> <li>Human Immunodeficiency virus [HIV-1] [AIDS]</li> <li>Influenza virus type A2 [Flu]</li> <li>Norovirus</li> <li>Poliovirus serotype 1 [PV-1]</li> <li>Respiratory Syncytial virus [RSV]</li> <li>Rhinovirus</li> <li>Rotavirus Strain WA</li> </ul>	<ul style="list-style-type: none"> <li>Mycobacterium bovis</li> </ul>	<ul style="list-style-type: none"> <li>Trichophyton mentagrophytes (Athletes Foot Fungus)</li> </ul>



airtechniques.com

## Order Information

Order No.	Product	Description
H6171	Monarch Surface Disinfectant Wipes	100 – 7"x9" wipes per canister (12 canisters per shipper case)
H6186	Monarch Surface Disinfectant Wipes	160 – 6"x6.75" wipes per canister (12 canisters per shipper case)
H6110-N	Monarch Surface Disinfectant Spray	32 oz. spray bottle (12 bottles per shipper case)



CORPORATE HEADQUARTERS  
1295 Walt Whitman Road  
Melville, New York 11747- 3062  
800-247-8324

Monarch is a registered trademark  
of Air Techniques, Inc.  
© Copyright 2018  
PN 9922-28 Rev. D

**AIR**  
**TECHNIQUES** equipped for life®