

The new, sharper claris *i4D*

The ultimate everyday intraoral camera.
Now clearer images, LED on/off feature,
and secure USB connection.



claris i4D

DIGITAL INTRAORAL CAMERA

CLARITY THROUGH TECHNOLOGY

BRILLIANT

- Highest quality optics for crystal clear images
- 8 white LEDs

ERGONOMIC

- Ergonomic handpiece that weighs only 3.4 oz
- Integrates seamlessly with imaging software

PRECISE

- Precision controlled sliding focus
- Silicone capture button for image stability
- Light on/off button function

SECURE

- Secure locking USB connection
- High grade shielded USB 2.0 cable

TABLET PC READY



Try our **TOUCH SCREEN POWERED**
cliosoft imaging software

*The tablet PC and ClioSoft software are sold separately.


Specifications

Imager	Sony Super HAD Color 1/4" CCD
Resolution	720 X 480 Pixels
Frame Rate	30 Fps Max.
Software Interface	Twain, WDM
Data Format	YUV 4:2:2, ITU-R 656
Field of View	90 Degrees
Focus Range	4mm-infinity
Focus Type	Variable Focus
Lgith Source	8 White LEDs
Signal to Noise Ratio	51 Db
Video Output	USB 2.0 Digital
Power Requirements	5v Dc, Powered By USB Port
Weight	2.4 Oz
Dimensions	8.3(L) X 0.85(W) X 1.15(H) Inches

System Requirements

- Intel Pentium 4 1.5GHz or higher
- 256MB of RAM or higher
- XVGA (1024 x 768) or higher resolution
- Video adapter and monitor
- CD-ROM or DVD drive
- Available USB 2.0 port
- Microsoft Window 2000/XP or later

Works with 32-bit & 64-bit versions of



T: 1. 888. 266. SOTA
F: 1. 714. 532. 6107
1073 N. Batavia St. Suite B
Orange, CA 92867

www.sotaimaging.com