

Digital Water Bath

With Variable Temperature Control
For Dental Operatory Use



A thermostatically controlled unit that maintains temperatures in a range of 32° C to 71° C (90° F to 160° F) to soften a variety of wax or compound materials.

- Digital temperature control with indicator light
- Removable stainless steel pan for easy sterilization
- Compatible with Hanau Water Bath pans
- LED screen which constantly displays the temperature in °F or °C
- Temperature indicator light shows when unit is heating
- Safety feature: unit stops heating when pan is removed
- UL and CUL approved

TECHNICAL SPECIFICATIONS

Voltage	AC 115/230V 50/60 Hz
Dimensions	19 cm (7.5") Wide 12.7 cm (5.0") High 22.2 cm (8.75") Deep
Capacity (Pan)	16.5 cm (6.5") Wide 10.8 cm (4.25") High 7.0 cm (2.75") Deep
Temperature	32° C–71° C (90° F–160° F) Variable
Weight	2.3 kg (5 lbs.)

ITEM NO. DESCRIPTION

05350 (115V)	Water Bath including removable stainless steel pan. (Lift Tray available separately)
05335 (230V)	Water Bath including removable stainless steel pan. (Lift Tray available separately)

ACCESSORIES

05327	#1576A Removable Pan
05326	#1578 Lift Tray



more reasons. one source

Whip Mix Corporation • 361 Farmington Ave. • Louisville, KY 40209
502-637-1451 • 800-626-5651 • www.whipmix.com