



# The Light & The Light 405


CURING LIGHTS



## EXPERIENCE NEW LIGHT TECHNOLOGY



### Comparison with Other Lights




## TOUCH of GENIUS™

- Concentrated light provides deep and efficient polymerization
- Reduce eye fatigue by removing unnecessary dispersed light
- Easy to read OLED screen
- Variety of curing modes
- Lightweight handpiece


Dimensions (Includes Battery)	232 x 26 x 26 cm
Weight	130 g
Wave length (The Light)	440~490 nm
Wave length (The Light 405)	405~490 nm
Power source	Lithium-ion Battery
Light source	High power LED

## Two Lights - One Great Result Every Cure





- Evenly dispersed concentrated light
- Time setting available
- Over 400 cures in one charge (high mode)
- Over 1,000 cures in one charge (low mode)
- Interchangeable battery
- 330° head twist



**Competitor  
LED - insufficient  
light dispersion.**



**The Light  
- uniform  
concentrated  
light increases  
polymerization  
performance.**



**The Light 405 has a dual wavelength used to cure all light cured restorative materials, as well as OPTIGLAZE™ color, our light-cured characterization coating for composite and acrylic indirect restorations.**

With OPTIGLAZE™ color, you can now change shades, shift value, characterize, create effects and glaze all your resin type prosthetics in no time at all. Best of all, offer your dentists and patients a more beautiful and esthetic restoration that will last due to GC's innovative new homogeneously dispersed nano-filler technology.

**402400 The Light (Operator)**  
Contains: One hand piece, one charger, one battery, one user manual and one box of 100 disposable covers.

**402402 Replaceable Battery**  
Contains: One battery.

**402405 The Light (Operator & Laboratory)**  
Contains: One hand piece, one charger, one battery, one user manual and one box of 100 disposable covers.

**402403 Disposable Covers**  
Contains: 100 replaceable disposable covers.

**008408 OPTIGLAZE™ color Set**  
Contains: All color shades (2.6 mL each), disposable palette, flat and round brush and brush holder.

