

NU SOURCE

INTERNATIONAL • INC.

PRODUCT FEATURES & BENEFITS

ABX PRE-MIXED Developer & Fixer

Packaged, 2 Gal. Developer & 2 Gal. Fixer

Specifically designed for all types and brands of dental x-ray film, as well as all brands of automatic film processors.

PERI-DENT PRE-MIXED Developer & Fixer

Packaged, 3 Qts. Developer & 3 Qts. Fixer

Specifically designed for all types of rollerless film processors and manual processing tanks.

FEATURES

BENEFITS

COLOR - ✓ Color Coding System

Specially formulated to help eliminate possible cross contamination. If cross contamination occurs when using the Red developer and Blue fixer, the contaminated product becomes clear.

FIXSCENT , No Odor Formulation

A specially formulated deodorizer completely eliminates chemical odors and provides a fresh pleasant floral fragrance.

Excellent Resistance to Oxidation

New formulation that resist oxidation (works cleaner and has longer active life).

High Capacity

New formula provides longer active life of product. It will develop more films over a longer period while maintaining high quality images.

INSTA-DENT Rapid X-Ray Developer & Fixer

Packaged, 2 Gal. Developer & 2 Gal. Fixer or
3 Qts Developer & 3 Qts Fixer.

Ready to use Developer & Fixer for instant processing of all X-Ray films. For chair side Endo processing and rapid manual hand tank processing.

FEATURES

BENEFITS

FIXSCENT, No Odor Formulation

A specially formulated deodorizer completely eliminates chemical odors and provides a fresh pleasant floral fragrance.

LOW OXIDATION Formula

Unique formulation resist oxidation and will not turn black in cups or tanks. Clean working to last twice as long and provide excellent image quality.

LH, LOW HAZARD SYSTEMS CLEANER

Packaged, 6 Qts. per case.

A one part **FAST ACTING** concentrate powder systems cleaner, that is designed to speedily dissolve silver oxide deposits and other stains from developer and fixer tanks and roller assemblies in automatic processors, and manual tanks.

FEATURES

BENEFITS

LOW HAZARD

Specially formulated to be personally and environmentally safe. No toxic fumes or odors.

CONCENTRATE FORMULA

A one part concentrate, just add water and it's ready to use. Works on both developer and fixer rollers or hand tanks. Mix full strength in quart container or dilute up to two gallons for soaking.

FAST ACTING

Designed to work on contact to dissolve silver oxide deposits and other stains from developer and fixer roller assemblies in automatic processors and manual tanks. If soaking is required, soak for no more than 20 to 30 minutes.

NO HAZARDOUS DUTY FREE

LH Systems Cleaner is environmentally safe and requires NO HAZARDOUS DUTY FEE when shipping or receiving.

MDX Manual Dental X-ray Concentrate Developer & Fixer

Packaged, 3 Qts. Developer/Fixer per case

Easy to use one part liquid concentrate for manual processing of all dental/medical x-ray films.

FEATURES

BENEFITS

FIXSCENT, No Odor Formulation

A specially formulated deodorizer completely eliminates chemical odors and provides a fresh pleasant floral fragrance.

Excellent Resistance to Oxidation

New formulation that resist oxidation (works cleaner and has longer active life).

High Capacity

New formula provides longer active life of product. It will develop more films over a longer period while maintaining high quality images.

For further information on these and other fine products, contact

NU SOURCE
INTERNATIONAL • INC.

WHITE MOUNTAIN IMAGING
46 CHENELL DRIVE, CONCORD NH 03301
(603)648-2124