

We Have Amalgam Separation Down to a Science

Solmetex has quickly become the industry leader in Amalgam Separation. What makes Solmetex different is our “Total Solution Approach” — integrating the simple design of the NXT Hg5 with waste handling and recycling into our complete product line, providing a truly green set of solutions for the dental practice.



Solutions

The NXT Hg5 provides the complete solution:

- New compact design, internal manifold, ideal for hard to fit areas
- Functional with all wet and dry vacuums
- Environmentally friendly packaging
 - Simple mail-back recycling program
 - Meets all “EPA Dental Rule” requirements
 - ISO 11143:2008 certified by NSF International



Easy

- Clear system, easy to inspect
- Two-minute container change out
- No tools necessary

Flexible

- Versatile, wall or floor mount
- Compatible with wet or dry vacuums
- Three available sizes to fit your practice

Support

- Compliance center
- Online certificates 24/7
- AutoShip program

Recycle

- Ship container via UPS
- EPA certified mercury facility
- Recycle Kit includes container and recycle packaging

Ordering Information

To order your Solmetex NXT Hg5 System, please contact your local dental product dealer. If you do not have a dealer, please visit www.solmetex.com and click on “dealer locator” to find a dental product dealer near you or contact Solmetex for information at **1-800-216-5505**.

NXT Hg5 Series of Amalgam Separators



**NXT Hg5 mini
Amalgam Separator**

Item # **NXT-HG5-MINI**

Specifications:

- 25”h x 11”w x 8”d
- ¾” inlet & outlet piping
- Fernco® couplers
- Floor / Wall mount bracket
- Services 1 - 4 chairs
- Container Reorder No. **NXT-HG5-002CR**



**NXT Hg5
Amalgam Separator**

Item # **NXT-HG5-001**

Specifications:

- 30”h x 11”w x 8”d
- 1½” inlet & outlet piping
- Fernco® couplers
- Floor / Wall mount bracket
- Services 1 - 10 chairs
- Container Reorder No. **NXT-HG5-002CR**



NEW

**NXT Hg5 High Volume
Amalgam Separator**

Item # **NXT-HG5-HV**

Specifications:

- 29¾”h x 17⅝”w x 22½”d
- 3” inlet & outlet piping
- Wall mount bracket
- Services 11 - 20 chairs
- Container Reorder No. **NXT-HG5-002CR**

Hg5 Series of Collection Containers



**NXT Hg5
Collection Container**

with Recycle Kit

Item # **NXT-HG5-002CR**



**Hg5 [Legacy]
Collection Container**

with Recycle Kit

Item # **HG5-002CR**

All Collection Container with Recycle Kits contain:

- Collection Container
- Proper packaging for shipment to our EPA approved recycling facility
- Recycling documentation available online 24/7
- *No contracts!*

Hg5 Series of Compliance Kits



**NXT Hg5
Compliance Kit**

Item # **NXT-HG5-CK**



**Hg5 [Legacy]
Compliance Kit**

Item # **HG5-CK**

All Compliance Kits contain:

- Collection Container
- Practice Waste Solutions Amalgam Bucket
- Proper packaging for shipment to our EPA approved recycling facility
- Recycling documentation available online 24/7

Solmetex -- Clearly the leader in amalgam separation. The NXT Hg5 amalgam separator is the industry standard and has won more awards than all competitors combined.

