



REXTAR X

PORTABLE XRAY GENERATOR



Registration Number K132041



www.gooddrs.us

1-844-448-5050

info@gooddrs.us

17609 Ventura Blvd. #110, Encino CA 91316

Lightweight, Safe & Reliable

Portable Handheld X-Ray Generator

1 Features of Rextar-X



Portable Handheld

High speed digital handheld, extremely mobile. One Rextar X can be used in multiple rooms. Saves on un-needed equipment cost.



Smart Buttons

Select Adult or Child patient, tooth type, receptor type (Film, Phosphor plates, or Digital Sensors) to provide a recommended dosage exposure time. Memorizes your favorite settings.



Safe to Use

Rextar X technology provides the highest level of safety. Double lead shielding and scatter shield protects the operator and patient. Nearly half the exposure time of other handhelds.



Compact & Light

Just over 5lbs, smaller and lighter than most existing portable X-Ray handheld devices. Camera-like design with neck strap.



Large Battery Capacity

Up to 300 exposures on one battery charge.


2 Specifications

XRy Generator	Type	High Frequency Inverter
	Power Output	140W
	Power Requirement	11.1 VDC
	Weight	Approx. 5.5 lbs
	Tube mA / kV	70kV/ 2mA
XRy Tube	Battery Charger	100-240 VAC, 50-60Hz 1.7A
	Model Name	Toshiba D- 041 Oil Cooled
	Inherent Filtration	>1.5 Al equivalent
XRy Control	Focal Spot	0.4 mm
	Exposure Time	Stationary

6 Compatible with All Sensors

- The Rextar X is compatible with all forms of X-ray detectors (Conventional Film, Phosphor Plates, and Digital Sensors)

6 Sample Images



Largest LCD Display Screen