

OPTIM 33TB, exceptional cleaning & disinfection capabilities.

A breath of fresh air.



- One step disinfectant cleaner – fast & convenient
- AHP – patented technology designed to deliver better cleaning and faster disinfection results (active ingredient 0.5% hydrogen peroxide)
- Hydrogen peroxide does not leave behind active residue which may lead to the development of resistant strains of microorganisms



Your Infection Control Specialist™

SciCan
A Sanavis Group Company



Corporate Member



FAST CONTACT TIME:

- 1 minute bactericidal & virucidal claim
- 30 second sanitizer
- 5 minute TB claim



- No strong chemical odor
- Active ingredient breaks down into oxygen and water
- Non-flammable
- EPA Cat. 4 for skin & eye irritation
- 0 HMIS health hazard rating*
- EPA Cat. 4: lowest toxicity level on EPA toxicity category

Active Ingredient	EPA Toxicity Category	Signal Word Required
AHP	Cat.4	NONE
Quats+Low % Alcohol	Cat.3	CAUTION
Quats+High % Alcohol	Cat.2	WARNING
Phenolics	Cat.1	DANGER

* Refer to HMIS rating on MSDS & compare with our competitors' wipes

SURFACE DISINFECTANT CONTACT TIME**

Product Name	Active Ingredient	Mycobacteria		Non-enveloped virus				Enveloped virus				Bacteria			
		TB	Poliovirus	Norovirus	Rotavirus	HIV	HBV (Hep.B)	HCV (Hep.C)	Influenza A2	MRSA	VRE	E.coli	Pseudomonas		
OPTIM 33 TB	AHP	5	1	1	1	1	1	1	1	1	1	1	1		
Birex SE	Phenolics	10	10			10	10		10	10		10	10		
Prospray	Phenolics	10	10			1	10	10	10	10		10	10		
Caviwipes	Quats-alcohol	3				2	2	2		2	2		3		
Opticide 3	Quats-alcohol	3	3			3	3	3	3	3	3	3	3		
Sanicloth Plus	Quats-alcohol	3				2	2	2		2	2	2	3		
Lysol IC	Quats-alcohol	10	10			10							10		

** Contact time in minutes. Contact time was obtained from publicly available material & is valid when table was constructed. Information will be updated if there are changes in claims. If blank, no published information was available.

ORDERING INFORMATION FOR OPTIM 33TB	
Ready-To-Use Pre-saturated Wipes:	Liquid:
OPT33-W12 RTU Wipes 6" x 7" (12 x 160 sheets/case)	OPT33-12x32 1 case of 12 x 32 oz.
Dry Wipes & Dry Wipes Refill:	OPT33-4x1 1 case of 4 x 1 gallon
OPT-EMPCAN STAT-Wipes Canister & Roll of Dry Wipes (12 x 160 sheets /case)	FSA01 Sprayer for 32 oz. bottle
SCI33R18 Dry Wipes Refill, 6" x 7" (18 x 160 sheets/case)	
SCI33R18XL Large Dry Wipes Refill 10" x 10" (18 x 60 sheets/case)	
OPT33-W10X10 RTU Wipes 10" x 10" (12 x 60 sheets/case)	

SciCan Inc., 701 Technology Drive, Canonsburg, PA 15317 • Phone (724) 820-1600; Fax (724) 820-1479; Toll Free 1-800-572-1211
www.scican.com



Your Infection Control Specialist™

OPTIM is a registered trademark and "Your Infection Control Specialist" is a trademark of SciCan Ltd. AHP is a registered trademark of Virox Technologies Inc. used under license. SD-341 US EN R3