

SAFETY DATA SHEET

Issuing Date 28-Aug-2019 Revision date 20-Dec-2019 Revision Number 0

1. Identification

Product identifier

Product Name Sani-Cloth HB Germicidal Disposable Wipe

Other means of identification

Product Code(s) SDS 0042-00 ENGLISH

Synonyms None

Registration Number(s) EPA Reg. No. 61178-4-9480

Other information 4BQ08401

Recommended use of the chemical and restrictions on use

Recommended useUse as a disinfectant on hard, non-porous surfaces. Read and understand the entire label

before using. Use only according to label directions. It is a violation of Federal law to use

this product in a manner inconsistent to label directions.

Restrictions on use For professional and hospital use

Details of the supplier of the safety data sheet

Manufacturer Address

Professional Disposables International, Inc. 400 Chestnut Ridge Road

Woodcliff Lake, NJ 07677

Emergency telephone number

Emergency Telephone PERS:

1-800-633-8253 (Domestic/Canada) 1-801-699-0667 (International)

2. Hazard(s) identification

Classification

Appearance Liquid saturated on wipe Physical state Liquid Odor Neutral/Nondescript

Label elements

Hazard statements

Not classified.

Precautionary Statements - Prevention

Observe good industrial hygiene practices

Precautionary Statements - Response

Wash hands after handling.

Precautionary Statements - Storage

Store in accordance to package label instructions.

Store away from incompatible materials.

Precautionary Statements - Disposal

Dispose of waste and residues in accordance with local authority requirements.

Other information

None known.

3. Composition/information on ingredients

Substance

Not applicable.

Mixture

Chemical name	CAS No.	Weight-%	Hazardous Material Information Review Act registry number (HMIRA registry #)	Date HMIRA filed and date exemption granted (if applicable)
Tetrasodium EDTA	64-02-8	0.1384	-	-

4. First-aid measures

Description of first aid measures

General advice If you feel unwell, seek medical advice (show the label where possible). In case of a severe

allergic or anaphylactic reaction, contact Poison Control immediately.

Inhalation Remove to fresh air.

Eye contact Immediately flush with plenty of water. After initial flushing, remove any contact lenses and

continue flushing for at least 15 minutes.

Skin contact IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin

with water/shower. (Call a physician if symptoms occur).

Ingestion Do not induce vomiting unless directed to do so by a medical professional. Never give

anything by mouth if victim is unconscious or is convulsing.

Most important symptoms and effects, both acute and delayed

Symptoms No information available.

Indication of any immediate medical attention and special treatment needed

5. Fire-fighting measures

Suitable Extinguishing Media Use extinguishing measures that are appropriate to local circumstances and the

surrounding environment.

Unsuitable extinguishing media CAUTION: Use of water spray when fighting fire may be inefficient.

Specific hazards arising from the

chemical

Firefighters should wear a self-contained breathing apparatus.

Explosion data

Sensitivity to mechanical impact None.

Sensitivity to static discharge None.

Special protective equipment for

fire-fighters

Firefighters should wear self-contained breathing apparatus and full protective clothing for fires involving chemicals or in confined spaces. Use personal protection equipment.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

Do not touch damaged containers or spilled material unless wearing appropriate protective Personal precautions

clothing. Follow standard industrial hygiene practices. Wear appropriate protective

equipment and clothing during clean-up. For personal protection, see section 8 of the SDS.

Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so. Stop leak if you can do it without risk.

Cover damaged package or spilled material with damp towel or rag and keep wet with liquid Methods for cleaning up

bleach or other disinfectant.

Prevention of secondary hazards Environmental precautions: Avoid discharge into lakes, streams, ponds or public waters.

See section 13 for more information. Reference to other sections

7. Handling and storage

Precautions for safe handling

Handle in accordance with good industrial hygiene and safety practice. Advice on safe handling

Conditions for safe storage, including any incompatibilities

Storage Conditions Store in a cool, dry location away from incompatible materials. Do not store near heat or

> open flame. Do not freeze or expose to extreme cold for a long period. Do not contaminate water, food or feed by storage or disposal. For containers: Protect container from physical

damage. When not in use keep center cap of lid closed to prevent moisture loss.

8. Exposure controls/personal protection

Control parameters

This product, as supplied, does not contain any hazardous materials with occupational **Exposure Limits**

exposure limits established by the region specific regulatory bodies.

Appropriate engineering controls

Engineering controls Showers

> Eyewash stations Ventilation systems.

Individual protection measures, such as personal protective equipment

Eve/face protection No special protective equipment required under normal use conditions. If needed, wear eye

protection per facility protocol to avoid eye contact.

Skin and body protection No special protective equipment required under normal use conditions. If needed, wear

protective gloves per facility protocol to avoid skin contact.

No protective equipment is needed under normal use conditions. If exposure limits are Respiratory protection

exceeded or irritation is experienced, ventilation and evacuation may be required.

Handle in accordance with good industrial hygiene and safety practice. General hygiene considerations

9. Physical and chemical properties

Information on basic physical and chemical properties

Physical state Liquid

Appearance Liquid saturated on wipe

Color White

Odor Neutral/Nondescript **Odor threshold** No information available

Remarks • Method Property Values

9.8-11.3 (Liquid) pН

No data available Melting point / freezing point None known Boiling point / boiling range No data available None known Flash point Tag Closed Cup >100.00 °C / >212 °F **Evaporation rate** No data available None known Flammability (solid, gas) No data available None known Flammability Limit in Air None known

Upper flammability or explosive No data available

limits

Lower flammability or explosive No data available

limits

No data available None known Vapor pressure Vapor density No data available None known Relative density No data available None known No data available Water solubility None known Solubility in other solvents Wipe is not soluble None known Partition coefficient No data available None known **Autoignition temperature** No data available None known **Decomposition temperature** No data available None known Kinematic viscosity No data available None known **Dynamic viscosity** No data available None known

Other information

No information available **Explosive properties Oxidizing properties** No information available Softening point No information available Molecular weight No information available

Disposable Wipe

Revision date 20-Dec-2019

VOC Content (%)

Liquid Density

Bulk density

No information available
No information available

10. Stability and reactivity

Reactivity May react with incompatible materials.

Chemical stability Stable under normal conditions.

Possibility of hazardous reactions Hazardous polymerization does not occur.

Conditions to avoid Avoid temperatures exceeding the flash point.

Incompatible materials Acids. Oxidizing agent. Caustics.

Hazardous decomposition products May include and are not limited to: Oxides of carbon. Oxides of phosphorus.

11. Toxicological information

Information on likely routes of exposure

Product Information

Inhalation Not a normal route of harmful exposure.

Eye contact Specific test data for the substance or mixture is not available.

Skin contact Specific test data for the substance or mixture is not available.

Ingestion Specific test data for the substance or mixture is not available.

Symptoms related to the physical, chemical and toxicological characteristics

Symptoms No information available.

Acute toxicity

Numerical measures of toxicity

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
Tetrasodium EDTA 64-02-8	= 10 g/kg (Rat) = 1658 mg/kg (Rat)	-	-

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation Prolonged skin contact may cause temporary irritation.

Serious eye damage/eye irritation Direct contact with eyes may cause temporary irritation.

Respiratory or skin sensitization This product is not expected to cause skin sensitization.

Germ cell mutagenicity No information available.

Carcinogenicity No information available.

Disposable Wipe

Reproductive toxicity No information available.

STOT - single exposure No information available.

STOT - repeated exposureNo information available.

Aspiration hazard No information available.

12. Ecological information

Ecotoxicity

Chemical name	Algae/aquatic plants	Fish	Toxicity to	Crustacea
			microorganisms	
Tetrasodium EDTA	EC50: =1.01mg/L (72h,	LC50: =41mg/L (96h,	-	EC50: =610mg/L (24h,
64-02-8	Desmodesmus	Lepomis macrochirus)		Daphnia magna)
	subspicatus)	LC50: =59.8mg/L (96h,		
		Pimephales promelas)		

Persistence and degradability No information available.

Bioaccumulation No information available.

Other adverse effects No other adverse environmental effects (e.g. ozone depletion, photochemical ozone

creation potential, endocrine disruption, global warming potential) are expected from this

Revision date 20-Dec-2019

component.

13. Disposal considerations

Waste treatment methods

Waste from residues/unused

products

Dispose of in accordance with local regulations. Dispose of waste in accordance with

environmental legislation. Do not flush in toilet.

Contaminated packaging Do not reuse empty containers.

14. Transport information

DOT Not regulated as dangerous goods.

TDG Not regulated as dangerous goods.

MEX Not regulated as dangerous goods.

ICAO (air) Not regulated as dangerous goods.

IATA Not regulated as dangerous goods.

IMDG Not regulated as dangerous goods.

RID Not regulated as dangerous goods.

ADR Not regulated as dangerous goods.

ADN Not regulated as dangerous goods.

15. Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture

International Regulations

The Montreal Protocol on Substances that Deplete the Ozone Layer Not applicable

The Stockholm Convention on Persistent Organic Pollutants Not applicable

The Rotterdam Convention Not applicable

International Inventories

Contact supplier for inventory compliance status. **TSCA**

DSL/NDSL Does not comply **EINECS/ELINCS** Does not comply

Contact supplier for inventory compliance status. **ENCS IECSC** Contact supplier for inventory compliance status. **KECL** Contact supplier for inventory compliance status. **PICCS** Contact supplier for inventory compliance status. Contact supplier for inventory compliance status. **AICS**

Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

US Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

SARA 311/312 Hazard Categories

Should this product meet EPCRA 311/312 Tier reporting criteria at 40 CFR 370, refer to Section 2 of this SDS for appropriate classifications. Under the amended regulations at 40 CFR 370, EPCRA 311/312 Tier II reporting for the 2017 calendar year will need to be consistent with updated hazard classifications.

CWA (Clean Water Act)

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

CERCLA

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material.

US State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

U.S. State Right-to-Know Regulations

U.S. EPA Label Information KEEP OUT OF REACH OF CHILDREN CAUTION Precautionary statements Hazards to humans and

domestic animals Caution. Causes moderate eye irritation. Harmful if absorbed through skin. Avoid contact with eyes, skin or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum or tobacco. Remove contaminated clothing and wash before reuse.

EPA Pesticide Registration Number 61178-4-9480

16. Other information

NFPA Health hazards 0 Flammability 0 Instability 0 Physical and chemical

properties -

HMIS Health hazards 0 Flammability 0 Physical hazards 0 Personal protection X

Key or legend to abbreviations and acronyms used in the safety data sheet

Legend Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

TWA TWA (time-weighted average) STEL STEL (Short Term Exposure Limit)

Ceiling Maximum limit value Skin designation

Key literature references and sources for data used to compile the SDS

Agency for Toxic Substances and Disease Registry (ATSDR)

U.S. Environmental Protection Agency ChemView Database

European Food Safety Authority (EFSA)

EPA (Environmental Protection Agency)

Acute Exposure Guideline Level(s) (AEGL(s))

U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act

U.S. Environmental Protection Agency High Production Volume Chemicals

Food Research Journal

Hazardous Substance Database

International Uniform Chemical Information Database (IUCLID)

Japan GHS Classification

Australia National Industrial Chemicals Notification and Assessment Scheme (NICNAS)

NIOSH (National Institute for Occupational Safety and Health)

National Library of Medicine's ChemID Plus (NLM CIP)

National Library of Medicine's PubMed database (NLM PUBMED)

National Toxicology Program (NTP)

New Zealand's Chemical Classification and Information Database (CCID)

Organization for Economic Co-operation and Development Environment, Health, and Safety Publications

Organization for Economic Co-operation and Development High Production Volume Chemicals Program

Organization for Economic Co-operation and Development Screening Information Data Set

RTECS (Registry of Toxic Effects of Chemical Substances)

World Health Organization

Issuing Date 28-Aug-2019

Revision date 20-Dec-2019

Revision Note No information available.

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage,

transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

End of Safety Data Sheet