



MAKE A GREAT FIRST IMPRESSION WITH GC TRAYS

Impression
Trays



GC





TABLE OF CONTENTS

Disposable Trays.....	1-2
COE-Sideless Trays	1
CHECK-BITE™ Trays	1
COE Spacer Trays (Full-Arch Perforated)	2
COE Spacer Trays (Full-Arch Perforated Partial).....	2
Metal Trays.....	3-9
Nickel-Plated Check-Bite Trays	3
Nickel-Plated Check-Bite Trays (Complete Denture Perforated).....	3
Depressed Anterior (Perforated)	4
Extra Long (Perforated)	4
Mcgowan-Winkler Denture (Perforated).....	5
Mcgowan-Winkler Complete Denture (Perforated)	5
Pediatric (Perforated)	6
Regular (Perforated)	6
Partial (Perforated).....	7
Regular (Solid).....	7
Nickel-Plated Partial (Solid)	8
Stainless Steel (Perforated)	8
Stainless Steel (Perforated Partial)	9
Stainless Steel (Solid)	9
STO-K Metal Trays.....	10-11
Perforated	10
Solid (Aluminum)	11
COE TRAY PLASTIC	12
Adhesives	13
Tray Adhesive (COE ADHESIVE)	13
Tray Adhesive (COE ADHESIVE THINNER).....	13
Tray Adhesive (GC TRAY ADHESIVE REMOVER)	13
Tray Adhesive (UNIVERSAL VPS ADHESIVE)	13
Impression Materials	14
EXABITE™ II	14
EXACLEAR™	14
EXAFAST™	14
EXAFLEX®	14
EXAJET™	14
EXA'lence™	14
EXAMIX™	14
EQUIPMENT & ACCESSORIES	14
Notes.	15



DISPOSABLE TRAYS

COE-SIDELESS TRAYS

Sideless, two-part, bite registration trays are for use with all types of impression materials, can be used individually, or can be adjusted to four sizes to accommodate various size arches. Provides three impression functions in one easy step:

- Accurate impression
- Counter impression of opposing arch
- Precise bite registration

COE-SIDELESS

Triple-Function Impression Trays

Contains: Fifty trays per package.

150001 Sideless

Sizes		Width (cm)	Length (cm)	Depth (cm)
Tray	Sideless	2-1/2"	3"	N/A



COE-SIDELESS Triple-Function Impression Tray

CHECK-BITE™ TRAYS

These trays cover more than half the arch area, resulting in multiple bite registration points for more accurate articulation. Performs three impression techniques in one easy step:

- Accurate impression
- Counter impression of opposing arch
- Precise bite registration

CHECK-BITE™

Triple-Function Posterior Style

Contains: Fifty trays per box.

257350 #73D

Triple-Function Anterior Style

Contains: Forty trays per box.

257540 #75D

Sizes		Width (cm)	Length (cm)	Depth (cm)
73D	Posterior	2-3/8"	2-3/16"	3/4"
75D	Anterior	2-5/8"	1-1/2"	3/4"



Posterior #73D

Anterior #75D



DISPOSABLE TRAYS

COE SPACER™ TRAYS

A unique zig-zag spacer bar has been incorporated into the tray floor to minimize the amount of contact between the remaining teeth and the tray.

COE SPACER™ TRAYS - FULL ARCH (PERFORATED)

1 Shelf Pack - 72 Trays

Contains: Six bags (all one size).

- 250016 #1D (Large Upper)
- 250046 #4D (Medium Upper)
- 250076 #7D (Small Upper)
- 250206 #20D (Large Lower)
- 250216 #21D (Medium Lower)
- 250226 #22D (Small Lower)
- 250096 Assorted Package

1 Master Pack - 432 Trays

Contains: One carton of 6 shelf packs (all one size).

- 250012 #1D (Large Upper)
- 250412 #4D (Medium Upper)
- 250712 #7D (Small Upper)
- 252012 #20D (Large Lower)
- 252112 #21D (Medium Lower)
- 252212 #22D (Small Lower)
- 250112 Assorted Package

Upper Sizes		Width (cm)	Length (cm)	Depth (cm)
1D	Large	3"	2-1/2"	7/8"
4D	Medium	2-7/8"	2-3/8"	3/4"
7D	Small	2-5/8"	2"	11/16"
Lower Sizes		Width (cm)	Length (cm)	Depth (cm)
20D	Large	3"	2-3/8"	3/4"
21D	Medium	3"	2-1/4"	3/4"
22D	Small	2-15/16"	2"	13/16"



COE SPACER™ TRAYS - FULL ARCH (PERFORATED PARTIAL)

1 Shelf Pack - 72 Trays

Contains: Six bags (all one size).

- 250306 #30D (Upper Left, Lower Right)
- 250316 #31D (Upper Right, Lower Left)
- 250326 #32D (Anterior)

1 Master Pack - 432 Trays

Contains: One carton of 6 shelf packs (all one size).

- 253012 #30D (Upper Left, Lower Right)
- 253112 #31D (Upper Right, Lower Left)
- 253212 #32D (Anterior)

Sizes		Width (cm)	Length (cm)	Depth (cm)
30D	Upper Left or Lower Right	1-9/16"	2-1/4"	3/4"
31D	Upper Right or Lower Left	1-9/16"	2-1/4"	3/4"
32D	Anterior (Upper and Lower)	2-9/16"	1-3/8"	13/16"



METAL TRAYS

NICKEL-PLATED CHECK-BITE™ TRAYS

When poured and mounted on the articulator, CHECK-BITE™ metal double-arch impression trays produce a pair of complementary, separable models in their true occlusal relationship.

CHECK-BITE™ TRAYS

Unilateral (Water cooled - for Hydrocolloid)

267001 #70

Unilateral (Non-water cooled)

267201 #72

Occlusal Inserts

262025 Package of 25

Sizes		Width (cm)	Length (cm)	Depth (cm)
70	Water Cooled	3/4"	4-1/2"	1-1/16"
72	Non-water Cooled	3/4"	4"	15/16"



Unilateral
(Non-water cooled)
#72



Upper
#62



Lower
#67

COMPLETE DENTURE - PERFORATED

Note: The larger the tray number, the smaller the tray.

Set of 8

Contains Numbers: 61, 62, 63, 64, 66, 67, 68 and 69.

260608 Assorted Set

Individual Complete Denture Trays

260611 #61 (Upper)

260621 #62 (Upper)

260631 #63 (Upper)

260641 #64 (Upper)

260651 #66 (Lower)

260671 #67 (Lower)

260681 #68 (Lower)

260691 #69 (Lower)

		Width (cm)	Length (cm)	Depth (cm)
Upper Sizes				
61	Large	3-1/8"	2-3/4"	3/4"
62	Medium	3"	2-5/8"	5/8"
63	Small	2-3/4"	2-5/8"	9/16"
64	Smaller	2-11/16"	2-3/8"	1/2"
Lower Sizes				
66	Large	3-15/16"	2-1/2"	3/8"
67	Medium	2-3/4"	2-1/2"	3/8"
68	Small	2-7/8"	2-3/8"	3/8"
69	Smaller	2-1/8"	2-1/8"	5/8"



METAL TRAYS

DEPRESSED ANTERIOR - PERFORATED

Note: The larger the tray number, the smaller the tray.

Set of 6

Contains Numbers: A1, A4, A7, A20, A21 and A22.

268006 Assorted Set

Individual Depressed Anterior Trays - Perforated

268011 #A1 (Upper)

268041 #A4 (Upper)


268071 #A7 (Upper)

268201 #A20 (Lower)

268211 #A21 (Lower)

268221 #A22 (Lower)

Upper Sizes		Width (cm)	Length (cm)	Depth (cm)
A1	Large	3-3/16"	2-5/8"	3/4"
A4	Medium	2-13/16"	2-1/2"	3/4"
A7	Small	2-3/4"	2-1/4"	3/4"
Lower Sizes		Width (cm)	Length (cm)	Depth (cm)
A20	Large	3-3/16"	2-1/2"	3/4"
A21	Medium	3-1/8"	2-1/4"	3/4"
A22	Small	2-3/4"	2-5/16"	13/16"



EXTRA LONG - PERFORATED

Note: The larger the tray number, the smaller the tray.

Set of 8

Contains Numbers: XL5, XL7, XL9, XL14, XL21, XL22, XL24 and XL28.

269008 Assorted Set

Individual Complete Denture Trays

269051 #XL5 (Upper)

269071 #XL7 (Upper)

269091 #XL9 (Upper)

269141 #XL14 (Upper)

269211 #XL21 (Lower)

269221 #XL22 (Lower)

269241 #XL24 (Lower)

269281 #XL28 (Lower)

Upper Sizes		Width (cm)	Length (cm)	Depth (cm)
XL5	Large	2-3/4"	2-1/2"	3/4"
XL7	Medium	2-3/4"	2-3/8"	11/16"
XL9	Small	2-9/16"	2"	5/8"
XL14	Smaller	2-5/16"	1-7/16"	13/16"
Lower Sizes		Width (cm)	Length (cm)	Depth (cm)
XL21	Large	3-1/8"	2-5/8"	13/16"
XL22	Medium	2-3/4"	2-1/2"	13/16"
XL24	Small	2-11/16"	2-1/4"	11/16"
XL28	Smaller	2-1/8"	2-1/8"	5/8"

METAL TRAYS

McGOWAN-WINKLER DENTURE - PERFORATED

Note: The larger the tray number, the smaller the tray.

Set of 8

Contains Numbers: 41, 42, 43, 44, 45, 46, 47 and 48.

260408 Assorted Set

Individual McGowan - Immediate Denture Trays

- 260401 #40 (Upper)
- 260411 #41 (Upper)
- 260421 #42 (Upper)
- 260431 #43 (Upper)
- 260441 #44 (Upper)
- 260451 #45 (Lower)
- 260461 #46 (Lower)
- 260481 #48 (Lower)

Upper Sizes		Width (cm)	Length (cm)	Depth (cm)
40	Extra Large	3-5/16"	2-5/8"	7/8"
41	Large	3-1/8"	2-5/8"	7/8"
42	Medium	2-7/8"	2-7/16"	7/8"
43	Small	2-7/8"	2-1/2"	7/8"
44	Smaller	2-11/16"	2-3/8"	3/4"
Lower Sizes				
45	Large	2-7/8"	2-1/2"	7/8"
46	Medium	2-5/8"	2-3/8"	7/8"
48	Small	2-3/8"	2-7/16"	3/4"



McGOWAN-WINKLER COMPLETE DENTURE - PERFORATED

Note: The larger the tray number, the smaller the tray.

Set of 8

Contains Numbers: 80, 81, 82, 84, 85, 86, 87 and 88.

260808 Assorted Set

Individual McGowan-Winkler Complete Denture Trays

- 260801 #80 (Upper)
- 260811 #81 (Upper)
- 260821 #82 (Upper)
- 260841 #84 (Lower)
- 260851 #85 (Lower)
- 260861 #86 (Lower)
- 260871 #87 (Lower)
- 260881 #88 (Lower)
- 260891 #89 (Lower)

Upper Sizes		Width (cm)	Length (cm)	Depth (cm)
80	Large	3-1/8"	2-3/4"	3/4"
81	Medium	2-13/16"	2-1/2"	9/16"
82	Small	2-1/2"	2-3/8"	3/8"
Lower Sizes				
84	Extra Large	2-15/16"	2-3/4"	1/4"
85	Large	3-1/4"	2-7/16"	5/16"
86	Medium	3"	2-3/8"	1/4"
87	Small	3"	2-1/4"	1/4"
88	Smaller	2-3/4"	2-1/2"	3/8"
89	Extra Small	2-1/2"	2-5/16"	3/4"

Individual McGowan-Winkler Solid Trays

- 261871 #87



METAL TRAYS

PEDIATRIC - PERFORATED

Note: The larger the tray number, the smaller the tray.

Set of 6

Contains Numbers: 7, 9, 14, 22, 24 and 28.

260016 Assorted Set

Individual Pediatric Trays

260091 #9 (Upper)


260141 #14 (Upper)

260221 #22 (Lower)

260241 #24 (Lower)

260281 #28 (Lower)

Upper Sizes		Width (cm)	Length (cm)	Depth (cm)
9	Large	2-7/16"	1-7/8"	3/4"
14	Small	2-1/4"	1-5/8"	5/8"
Lower Sizes				
22	Large	2-15/16"	2-1/8"	13/16"
24	Medium	2-1/2"	1-15/16"	3/4"
28	Small	2-3/16"	1-3/4"	5/8"



REGULAR - PERFORATED

Note: The larger the tray number, the smaller the tray.

Set of 8

Contains Numbers: 1, 3, 4, 5, 7, 20, 21 and 22.

260018 Assorted Set

Individual Tray - Perforated

260011 #1 (Upper)

260031 #3 (Upper)

260041 #4 (Upper)

260051 #5 (Upper)

260071 #7 (Upper)

260911 #X1 (Upper)

260201 #20 (Lower)

260211 #21 (Lower)

260221 #22 (Lower)

260921 #X20 (Lower)

Upper Sizes		Width (cm)	Length (cm)	Depth (cm)
1	Extra Large	2-15/16"	2-3/4"	1/4"
3	Large	3-1/4"	2-7/16"	5/16"
4	Medium	3"	2-3/8"	1/4"
5	Small	3"	2-1/4"	1/4"
7	Smaller	2-3/4"	2-1/2"	3/8"
X1	Extra Large	3-3/8"	2-9/16"	7/8"
Lower Sizes				
20	Large	3"	2-3/8"	3/4"
21	Medium	3"	2-1/4"	7/8"
22	Small	2-15/16"	2-1/8"	13/16"
X20	Extra Small	3-3/8"	2-9/16"	7/8"

METAL TRAYS

PARTIAL - PERFORATED

Note: The larger the tray number, the smaller the tray.

Set of 4

Contains Numbers: 30, 31, 32 and 33.

260304 Assorted Set

Individual Trays (Partial Perforated)

260301 #30 (Upper Left, Lower Right)

260311 #31 (Upper Right, Lower Left)

260321 #32 (Anterior Lower)

260331 #33 (Anterior Upper)

260991 #99 (Swivel)

Sizes		Width (cm)	Length (cm)	Depth (cm)
30	Upper Left or Lower Right	1-3/4"	2-7/16"	3/4"
31	Upper Right or Lower Left	1-3/4"	2-7/16"	3/4"
32	Anterior Lower	2-1/2"	1-1/2"	13/16"
33	Anterior Upper	2-5/8"	1-7/16"	13/16"
99	Swivel	1-1/4"	2-3/16"	3/4"



REGULAR - SOLID

Note: The larger the tray number, the smaller the tray.

Set of 8

Contains Numbers: 101, 104, 107, 109, 120, 121, 122 and 124.

261018 Assorted Set

Individual Tray - Perforated

261011 #101 (Upper)

261031 #103 (Upper)

261041 #104 (Upper)

261051 #105 (Upper)

261071 #107 (Upper)

261911 #109 (Upper)

261201 #120 (Lower)

261211 #121 (Lower)

261221 #122 (Lower)

261241 #124 (Lower)

		Width (cm)	Length (cm)	Depth (cm)
Upper Sizes				
101	Extra Large	3"	2-1/2"	7/8"
103	Large	2-7/8"	2-1/2"	7/8"
104	Medium	2-7/8"	2-3/8"	7/8"
105	Small	2-5/8"	2"	3/4"
107	Smaller	2-5/8"	2"	11/16"
109	Extra Small	2-7/8"	1-7/8"	3/4"
Lower Sizes				
120	Large	3"	2-3/8"	3/4"
121	Medium	3"	2-1/4"	7/8"
122	Small	2-15/16"	2-1/8"	13/16"
124	Extra Small	2-1/2"	1-15/16"	3/4"



METAL TRAYS

NICKEL-PLATED PARTIAL (SOLID)

Note: The larger the tray number, the smaller the tray.

Set of 4

Contains Numbers: 130, 131, 132 and 133.

261304 Assorted Set

Individual Trays - Partial Solid

261301 #130 (Upper Left, Lower Right)

261311 #131 (Upper Right, Lower Left)

261321 #132 (Anterior Lower)

261331 #133 (Anterior Upper)

Swivel

261991 #199

Sizes		Width (cm)	Length (cm)	Depth (cm)
130	Upper Left or Lower Right	1-3/4"	2-7/16"	3/4"
131	Upper Right or Lower Left	1-3/4"	2-7/16"	3/4"
132	Anterior Lower	2-1/2"	1-1/2"	13/16"
133	Anterior Upper	2-5/8"	1-7/16"	13/16"
199	Swivel	1-1/4"	2-3/16"	3/4"



Swivel #199



STAINLESS STEEL - PERFORATED

Note: The larger the tray number, the smaller the tray.

Set of 8

Contains Numbers: S1, S3, S4, S5, S7, S20, S21 and S22.

264008 Assorted Set

Individual Trays - Perforated

264011 #S1 (Upper)

264031 #S3 (Upper)

264041 #S4 (Upper)

264051 #S5 (Upper)

264071 #S7 (Upper)

264911 #SX1 (Extra Large Upper)

264201 #S20 (Lower)

264211 #S21 (Lower)

264221 #S22 (Lower)

264921 #SX20 (Extra Large Lower)

		Width (cm)	Length (cm)	Depth (cm)
Upper Sizes				
S1	Large	3"	2-5/16"	13/16"
S3	Medium	2-3/16"	2-3/8"	13/16"
S4	Small	2-5/16"	2-1/8"	13/16"
S5	Smaller	2-3/4"	2-1/16"	3/4"
S7	Extra Small	2-11/16"	1-15/16"	3/4"
SX1	Extra Large	3-7/16"	2-1/2"	13/16"
Lower Sizes				
S20	Large	3"	2-3/8"	13/16"
S21	Medium	2-7/8"	2-1/4"	13/16"
S22	Small	2-7/8"	2-3/16"	3/4"
SX20	Extra Large	3-1/4"	2-1/2"	13/16"



STO-K™ METAL TRAYS

STAINLESS STEEL - PERFORATED PARTIALS

Note: The larger the tray number, the smaller the tray.

Set of 4

Contains Numbers: S30, S31, S32 and S33.

264004 Assorted Set

Individual Trays - Perforated

264301 #S30 (Right Unilateral)

264311 #S31 (Left Unilateral)

264321 #S32 (Anterior Lower)

264331 #S33 (Anterior Upper)

Sizes		Width (cm)	Length (cm)	Depth (cm)
S30	Right Unilateral	1-1/4"	2-7/16"	3/4"
S31	Left Unilateral	1-1/4"	2-7/16"	3/4"
S32	Anterior Lower	2-5/8"	1-9/16"	13/16"
S33	Anterior Upper	2-5/8"	1-7/16"	3/4"



Anterior Upper
#S33



Right Unilateral
#S30



Upper
#S104



Lower
#S121

STAINLESS STEEL - SOLID

Note: The larger the tray number, the smaller the tray.

Set of 8

Contains Numbers: S101, S103, S104, S105, S107, S120, S121 and S122.

264108 Assorted Set

Individual Trays - Solid

264101 #S101 (Upper)

264103 #S103 (Upper)

264104 #S104 (Upper)

264105 #S105 (Upper)

264107 #S107 (Upper)

264120 #S120 (Lower)

264121 #S121 (Lower)

264122 #S122 (Lower)

		Width (cm)	Length (cm)	Depth (cm)
Upper Sizes				
S101	Extra Large	3"	2-5/16"	13/16"
S103	Large	2-13/16"	2-3/8"	13/16"
S104	Medium	2-5/16"	2-1/8"	13/16"
S105	Small	2-3/4"	2-1/16"	3/4"
S107	Smaller	2-11/16"	1-15/16"	3/4"
Lower Sizes				
S120	Large	3"	2-3/8"	13/16"
S121	Medium	2-7/8"	2-1/4"	13/16"
S122	Small	2-7/8"	2-3/16"	3/4"



STO-K™ METAL TRAYS

STO-K™ TRAYS

Anatomically designed to minimize pressure on lips, enhance patient comfort and reduce post insertion adjustments. These basic shapes facilitate more accurate impressions.

PERFORATED

Set of 23

Contains: Ten Upper and 13 Lower.

263023 Assorted Set

Perforated - Square

263471 #U-1-S (Upper)

263472 #U-2-S (Upper)

263473 #U-3-S (Upper)

263481 #L-1-S (Lower)

263482 #L-2-S (Lower)

263483 #L-3-S (Lower)

263484 #L-4-S (Lower)

263485 #L-5-S (Lower)

Perforated - Tapered

263571 #U-1-T (Upper)

263572 #U-2-T (Upper)

263573 #U-3-T (Upper)

263581 #L-1-T (Lower)

263582 #L-2-T (Lower)

263583 #L-3-T (Lower)

263584 #L-4-T (Lower)

Perforated - Ovoid

263671 #U-1-O (Upper)

263672 #U-2-O (Upper)

263673 #U-3-O (Upper)

263674 #U-4-O (Upper)

263681 #L-1-O (Lower)

263682 #L-2-O (Lower)

263683 #L-3-O (Lower)

263684 #L-4-O (Lower)



Square
(Upper)
#U-2-S



Square
(Lower)
#L-2-S



Tapered
(Upper)
#U-2-T



Tapered
(Lower)
#L-2-T



Ovoid
(Upper)
#U-2-O



Ovoid
(Lower)
#L-2-O

Upper Sizes		Width (cm)	Length (cm)	Depth (cm)
U-1-S	Large	2-7/16"	1-3/4"	7/16"
U-2-S	Medium	2-1/2"	2"	3/4"
U-3-S	Small	3-1/4"	2-3/4"	7/8"
U-1-T	Large	2-1/2"	2-1/8"	9/16"
U-2-T	Medium	2-5/8"	2-1/4"	5/8"
U-3-T	Small	3-1/8"	2-11/16"	5/8"
U-1-O	Extra Large	2-1/2"	2"	9/16"
U-2-O	Large	2-5/8"	2-1/4"	5/8"
U-3-O	Medium	2-3/4"	2-3/8"	1/2"
U-4-O	Small	3"	2-1/2"	11/16"

Lower Sizes

		Width (cm)	Length (cm)	Depth (cm)
L-1-S	Extra Large	2-7/16"	1-3/4"	7/16"
L-2-S	Large	2-1/2"	2"	1/2"
L-3-S	Medium	2-7/16"	2-3/8"	1/2"
L-4-S	Small	2-11/16"	2-7/8"	3/8"
L-5-S	Smaller	3"	2-3/4"	3/8"
L-1-T	Extra Large	2-3/8"	2-9/16"	3/8"
L-2-T	Large	2-1/2"	2-5/16"	1/2"
L-3-T	Medium	2-7/16"	2-3/8"	1/2"
L-4-T	Small	3"	2-3/8"	3/8"
L-1-O	Extra Large	2-7/16"	2"	5/16"
L-2-O	Large	2-11/16"	2-1/8"	5/16"
L-3-O	Medium	2-5/16"	2-7/16"	5/16"
L-4-O	Small	2-15/16"	2-5/8"	1/2"

SOLID (ALUMINUM)

Set of 23

Contains: Ten Upper and 13 Lower.

261023 Assorted Set

Solid - Square

- 261471 #U-11-S (Upper)
- 261472 #U-12-S (Upper)
- 261473 #U-13-S (Upper)
- 261481 #L-11-S (Lower)
- 261482 #L-12-S (Lower)
- 261483 #L-13-S (Lower)
- 261484 #L-14-S (Lower)
- 261485 #L-15-S (Lower)

Solid - Tapered

- 261571 #U-11-T (Upper)
- 261572 #U-12-T (Upper)
- 261573 #U-13-T (Upper)
- 261581 #L-11-T (Lower)
- 261582 #L-12-T (Lower)
- 261583 #L-13-T (Lower)
- 261584 #L-14-T (Lower)

Solid - Ovoid

- 261671 #U-11-O (Upper)
- 261672 #U-12-O (Upper)
- 261673 #U-13-O (Upper)
- 261674 #U-14-O (Upper)
- 261681 #L-11-O (Lower)
- 261682 #L-12-O (Lower)
- 261683 #L-13-O (Lower)
- 261684 #L-14-O (Lower)



Upper Sizes		Width (cm)	Length (cm)	Depth (cm)
U-11-S	Large	2-7/16"	1-3/4"	7/16"
U-12-S	Medium	2-1/2"	2"	3/4"
U-13-S	Small	3-1/4"	2-3/4"	7/8"
U-11-T	Large	2-1/2"	2-1/8"	9/16"
U-12-T	Medium	2-5/8"	2-1/4"	5/8"
U-13-T	Small	3-1/8"	2-11/16"	5/8"
U-11-O	Extra Large	2-1/2"	2"	9/16"
	Large			
U-12-O	Large	2-5/8"	2-1/4"	5/8"
U-13-O	Medium	2-3/4"	2-3/8"	1/2"
U-14-O	Small	3"	2-1/2"	11/16"

Lower Sizes		Width (cm)	Length (cm)	Depth (cm)
L-11-S	Extra Large	2-7/16"	1-3/4"	7/16"
L-12-S	Large	2-1/2"	2"	1/2"
L-13-S	Medium	2-7/16"	2-3/8"	1/2"
L-14-S	Small	2-11/16"	2-7/8"	3/8"
L-15-S	Smaller	3"	2-3/4"	3/8"
L-11-T	Extra Large	2-3/8"	2-9/16"	3/8"
L-12-T	Large	2-1/2"	2-5/16"	1/2"
L-13-T	Medium	2-7/16"	2-3/8"	1/2"
L-14-T	Small	3"	2-3/8"	3/8"
L-11-O	Extra Large	2-7/16"	2"	5/16"
L-12-O	Large	2-11/16"	2-1/8"	5/16"
L-13-O	Medium	2-5/16"	2-7/16"	5/16"
L-14-O	Small	2-15/16"	2-5/8"	1/2"



COE TRAY PLASTIC™

COE TRAY PLASTIC™

Self-cured acrylic tray material for construction of custom impression trays. Available in Regular Set (15 minutes) and Fast Set (10 minutes).



COE TRAY PLASTIC™

Professional Package

Contains: Powder (1 lb; 454 g), Liquid (6 oz; 177 mL) and 1 COE-SEP™ (2 oz; 59 mL).

240111 Fast Set - White

240011 Regular Set - White

Economy Package

Contains: Powder (3 lbs; 1.4 kg), Liquid (1 pt; 473 mL) and 2 COE-SEP (4 oz; 118 mL).

240113 Fast Set - White

240013 Regular Set - White

Powder (Jar - 1 lb; 454 g)

240012 White

Powder (Jar - 3 lbs; 1.4 kg)

240014 White

Powder (Drum - 25 lbs; 11.4 kg)

240016 White

Liquid (Bottle - 6 oz; 117 mL)

240191 Fast Set

240091 Regular Set

Liquid (Can - 1 pt; 473 mL)

240192 Fast Set

240092 Regular Set

Liquid (Can - 1 qt; 946 mL)

240193 Fast Set

240093 Regular Set

Liquid (Gallon - 1 gal; 4-1 qt cans)

240194 Fast Set

240094 Regular Set



ADHESIVES

TRAY ADHESIVES

Regardless of the impression tray, use the tray adhesive made for the impression material (use for all types of trays). Without the adhesive, the impression:

- Shrinks toward the center of the mass
- Results in a smaller impression
- Crown may not seat completely



COE ADHESIVE*

133911 15 mL
133912 12 x 15 mL

COE ADHESIVE Thinner*

133921 15 mL

GC TRAY ADHESIVE REMOVER

250116 472 mL

UNIVERSAL VPS ADHESIVE

138991 7 mL
138912 12 x 7mL

* For use with COE-FLEX® and OMNIFLEX™ Impression Materials.



IMPRESSION MATERIALS

BITE REGISTRATION

EXABITE™ II

- 132412 Registration Creme (2-Pack with Tips; 48 mL)
- 132418 Registration Creme (8-Pack with Tips; 48 mL)

EXACLEAR™

- 008482 Translucent (2-Pack with Tips; 48 mL)

VPS

EXAFAST™

- 137007 Regular Body 2-Pack with Tips (48 mL)
- 137040 Regular Body Value 8-Pack with Tips (48 mL)
- 137107 Injection 2-Pack with Tips (48 mL)
- 137035 Injection Value 8-Pack with Tips (48 mL)
- 137207 Heavy Body 2-Pack with Tips (48 mL)
- 137275 Heavy Body 2-Pack with Tips (75 mL)
- 137295 Heavy Body Super Value 20-Pack (75 mL)
- 137030 Heavy Body Value 8-Pack with Tips (48 mL)
- 137407 Monophase 2-Pack with Tips (48 mL)
- 137045 Monophase Value 8-Pack with Tips (48 mL)
- 137305 PUTTY Clinic Package
- 137307 PUTTY Standard Package
- 137202 Intro Kit (2 Injection and 4 Heavy Body; 48 mL)

EXAMIX™

- 137016 Regular Body 2-Pack with Tips (48 mL)
- 137090 Regular Body Value 8-Pack with Tips (48 mL)
- 137106 Injection 2-Pack with Tips (48 mL)
- 137085 Injection Value 8-Pack with Tips (48 mL)
- 137206 Heavy Body 2-Pack with Tips (48 mL)
- 137080 Heavy Body Value 8-Pack with Tips (48 mL)
- 137406 Monophase 2-Pack with Tips (48 mL)
- 137095 Monophase Value 8-Pack with Tips (48 mL)

EXAFLEX®

- 138020 Regular Clinic Package
- 138901 Regular Standard Package
- 138101 Injection Standard Package
- 138120 Injection Clinic Package
- 138201 Heavy Standard Package
- 138220 Heavy Clinic Package
- 138401 Monophase Standard Package
- 138420 Monophase Clinic Package
- 138301 PUTTY Standard Package
- 138305 PUTTY Clinic Package

EXAJET™

- 800083 Heavy Body Fast Set Refill Package
- 800084 Heavy Body Fast Set Clinic Package
- 800085 Heavy Body Normal Intro Package
- 800086 Heavy Body Normal Refill Package
- 800087 Heavy Body Normal Clinic Package
- 800088 Cartridge (1 piece)
- 800089 Dynamic Mixing Tips (50)

VPES™

EXA'lence™

FAST SET

- 137902 Intro Kit
- 137014 Heavy Refill (48 mL)
- 137104 Heavy Rigid Refill (48 mL)
- 137204 Extra Light Refill (48 mL)
- 137304 Light Refill (48 mL)
- 137404 Monophase Refill (48 mL)
- 137175 Heavy Refill (75 mL)
- 137370 Heavy Standard Package (370 mL)
- 137470 Heavy Refill (370 mL)

REGULAR SET


- 137402 Intro Kit
- 137004 Heavy Refill (48 mL)
- 137224 Extra Light Refill (48 mL)
- 137334 Light Refill (48 mL)
- 137444 Monophase Refill (48 mL)
- 137375 Heavy Refill (75 mL)
- 137170 Monophase Standard Package (370 mL)
- 137270 Heavy Standard Package (370 mL)
- 137870 Heavy Refill (370 mL)
- 137970 Monophase Refill (370 mL)

PUTTY (VPS)

- 137405 Clinic Package
- 137805 Standard Package

EQUIPMENT & ACCESSORIES

- 439500 Renfert Sympress Mixing Machine
- 152101 GC CARTRIDGE DISPENSER II NDS (48 mL)
- 152102 GC CARTRIDGE DISPENSER II NDS (75 mL)
- 150100 Intraoral Tips Universal/NDS (100)
- 150300 Intraoral Mixing Tips (Yellow; 100)
- 153012 Mixing Tips Type 2L (Aqua; 12)
- 153048 Mixing Tips Heavy (Purple; 48)
- 153148 Mixing Tips Injection (Pink; 48)
- 153248 Mixing Tips (Aqua; 48)
- 159048 Mixing Tips (Yellow; 48)



EXA'lence™ Impressions - Photos courtesy of James Pavlatos, DDS

Innovative Impression Materials and Trays from GC

With EXA impression materials and GC trays, you get handling, accuracy, thixotropicity, elastic recovery, stackability and pourability at an exceptional level. The high hydrophilicity of EXA VPS materials means they flow uniformly, eliminate voids and capture details more precisely.

EXA impressional materials offer virtually no shrinkage (less than 0.2% in 24 hours), and are so dimensionally stable that an accurate model can be poured up to two weeks after the impression has been taken.

No matter what type of impression is being taken, EXA materials set quickly, flow well, and have the elastic recovery characteristics and tear strength that make the first impression a lasting impression.